

Analysis of cell protrusions using Fiji


Overview

- This protocol can be used to assess protrusion length in 2D fluorescence images of cells.
- Unlike traditional neurite tracing, this method does not require cells to be well separated, and does not require protrusions to be traceable to a cell body. The protrusion segments are detected and measured individually.
- The number of cell bodies in the image is also counted to assist with normalization to cell number.
- **Limitations:**
 - If protrusions cross, the length will be underestimated, because each segment up to the intersection is measured separately.
 - Cell bodies are included in the skeleton, leading to overestimation of protrusion length. However, they contribute a relatively small proportion of overall length in normal cells.

Setup

Fiji and custom scripts are needed for the image analysis.

1. Download and install **Fiji** from <https://imagej.net/software/fiji/>.
2. Download two custom scripts from <https://github.com/CUIMC-Confocal/zhiguo Zhang> :
 - 2.1. **analyze_protrusions.ijm**
 - 2.2. **prunebysize.bsh**
3. Move the script **prunebysize.bsh** into Fiji's *plugins/scripts* folder. To see this folder on Mac, right-click Fiji.app > *Show Package Contents*.




Data preparation


- Requirements/assumptions
 - Input images must be **2D** (not Z stack); they can be max projections of a Z stack.
 - Images must be **single timepoints**.
 - Images can be multichannel. You can specify which channel is used for the analysis.
 - Images can be any format that Fiji can read, including microscope formats such as .ND2.
 - The protocol was developed for images with a **pixel size of 0.37 μm** .
- Batch processing
 - The script will process all the images within a folder. Place all of the images you want to analyze in a single folder.

Image analysis

1. In Fiji, open the script ***analyze_protrusions.ijm***
2. In the Script Editor window, click **Run**.




3. Set the input and output folders, the file suffix, and the channel to analyze, then click OK.



4. When all images have been processed, the Fiji Log will show "Finished":



5. Check the designated output folder. It should contain 2 output images and 2 data tables for each input image, and a **batch summary**.



- 5.1. The output images are the **skeleton** (detected protrusions) and an **overlay** of the original fluorescence with the detected cell bodies and skeleton. These may be useful for visually validating the results.



- 5.2 The data tables for each image are the “skeleton info” and the “branch info.”

- The **branch info** contains the important length information about the individual protrusion segments. The Branch Length is used to calculate the results in the batch summary.
- The **skeleton info** relates to distinct skeletons in the image, and is not really informative if the cell processes overlap.

Skeleton ID	Branch length	V1 x	V1 y
1	88.296	9	76
1	69.983	63	139
1	27.556	113	200
1	16.899	132	236

	# Branches	# Junctions	# End-point voxels	# Junction voxels	# Slab voxels
1	9	4	2	8	166
2	23	11	6	17	678
3	4	2	1	3	20
4	50	24	19	52	1443
5	2	1	1	1	97

6. The **batch summary** contains the derived data for the whole folder (one row per input image).

- 6.1. **Total Length** is the sum of all skeleton segments in the image
- 6.2. **Number of cells** is the number of detected cell bodies
- 6.3. **Normalized length** is total length divided by the number of cells.
- 6.4. **Median branch length** is the median length of all skeleton segments in the image.




Image	Total Length	Number of cells	Normalized length	Median branch length
1 crop	3610.166	7	515.738	23.763
4 crop	2073.351	9	230.372	21.107

Technical details

- Cell counts are obtained by applying the default ImageJ threshold algorithm and applying binary Open and Watershed operations, then filtering by size.
- Protrusions are detected by a multi-step process:
 - [Local \(adaptive\) Phansalkar thresholding](#) to enhance dim processes
 - [Frangi vesselness](#) filter to enhance thin linear objects
 - Global [percentile thresholding](#) to detect the enhanced processes and cell bodies
 - Filtering to remove more circular objects
- The resulting objects, representing processes and their attached cell bodies, are skeletonized (thinned repeatedly until they are 1 pixel wide). This allows analysis of the protrusions as a network of segments (documentation: [Analyze Skeleton](#)).
 - The detected skeletons are "pruned" to remove the smallest segments (which are usually artifacts).
 - The pruned skeletons are analyzed to generate the values in the summary.

References: please cite these in your papers

- **Fiji:** Schindelin et al, "Fiji: an open-source platform for biological-image analysis," *Nature Methods* 2012, doi:10.1038/nmeth.2019.
- **Analyze Skeleton:** Arganda-Carreras, I., et al. (2010). 3D reconstruction of histological sections: Application to mammary gland tissue. *Microscopy Research and Technique*, 73(11), 1019–1029. doi:10.1002/jemt.20829
- **The CSMSR. Suggested text for acknowledgement:**
These studies used the Confocal and Specialized Microscopy Shared Resource of the Herbert Irving Comprehensive Cancer Center at Columbia University, funded in part through the NIH/NCI Cancer Center Support Grant P30CA013696.

Please contact us with any questions: Confocal_HICCC@cumc.columbia.edu

TS 11/11/2024